

Armfit+



Automatic Upper Arm Blood Pressure Monitor + ECG



Blood Pressure




ECG



Bluetooth

Track Your Health Easily

- **Blood Pressure Measurement:** systolic, diastolic blood pressure and pulse rate
- **ECG Measurement:** 30s recording with internal electrodes
- **Multiple Lead set:** Lead I, Lead II and Chest Lead
- **Sync data to App:** real-time monitoring and sharing report
- **Portable design:** lightweight cuff, OLED display and built-in rechargeable battery



Specification

Blood pressure measurement

Technology: Oscillometric Method
Pressure range: 0 - 300mmHg
Pressure measurement accuracy: ± 3 mmHg
Pulse rate range: 40 to 200 /min
Pulse rate accuracy: ± 2 /min

ECG recording

ECG length: 30s
Heart rate range: 30 to 250 bpm
Heart rate accuracy: ± 2 /min or $\pm 2\%$
whichever is greater

Physical

Size (main unit) : 135mm(L) \times 45mm(W) \times 20mm(H)
Weight (main unit): 240g
Cuff size: 22 - 42cm
Wireless connectivity: Bluetooth 4.0 BLE

Power supply

Battery type: Rechargeable lithium battery
Battery usage: 500 measurements

Storage

Blood pressure measurement : 50 records
ECG recording: 10 records



SHENZHEN CREATIVE INDUSTRY CO., LTD.

Website: en.lepumedical.com

Add: Floor 5, BLD 9, Balwangxin High-Tech Industrial Park, Songbai Road, Xili Street Nanshan District, 518110 Shenzhen, China

Email: info@creative-sz.com

Tel: + 86 755 2643 351 Fax: + 86 755 2643 0930