



# LF05A

## Oxygen Concentrator



Large-screen display

Background soundtrack available

Storage box, clean and convenient



Convenient rotary knob control

Atomizing interface &  
Oxygen outlet protection

Built-in water tank design



## Product Features

Imported molecular sieve from France, service life is 1.5 times longer than normal brand

First class air compressor, continuously work for 24hours\*7days

Condensate drainage Patent: Unique patent to extend the service life of molecular sieve

Multiple noise reduction patent: actual noise < 45 dB; quite when sleeping

0.1ml/min accurate control with voice broadcast

12 kinds of alarm function with remote control

## Technical Parameters

Oxygen concentration:	≥ 90% (V/V) (percent by volume at 5 L/min)
Oxygen discharge pressure:	30-60 kPa (maximum pressure at the oxygen outlet)
Atomizing rate:	0.12 ml/min - 0.90 ml/min
Laboratory-tested running noise:	45 dB(A)
Dimension:	450mm × 300mm × 550mm
Net weight:	≈ 16.5 KG
Rated voltage of oxygenator body:	a.c.220V 50Hz
Classification of electric shock degree:	Applied parts of BF type
Electromagnetic compatibility:	Class B, Group 1 according to GB4824
Classification of safety degree:	Non-AP/APG type. Not suitable when the flammable anesthetic gas is mixed with air, or oxygen or nitrous oxide
Model/specifications:	LF05A
Oxygen flow:	1-5 L/min (adjustable)
Input power:	450VA
Atomizing pressure:	30~120kPa
Capacity of storage tank:	About 2 ml-8 ml
Service life:	5 years
Electric shock protection classification:	Class II
Degree of protection against liquid inlet:	IPX0 unprotected equipment
Classification of running modes:	Continuous running
Software version:	V1.0

## Shenzhen lepu intelligent medical device co. LTD

📍 Add: Floor 3, BLD 9, BaiWangxin High-Tech Industrial Park,  
Songbai Road, Xili Street, Nanshan District, Shenzhen, China.

✉ E-mail: paul.han@lepu-medical.com

🌐 Web: www.lepucare.com

☎ Tel: 400-830-9392



# LBP70D

Sphygmomanometer





LEPU MEDICAL

LBP70D  
Sphygmomanometer



#### Cuff identification

Intelligent pressure, reducing measurement errors.



#### Live voice broadcast

Whole-course voice from measurement reminder to results analysis.



#### Lithium battery power supply

2200mAh lithium battery



#### 5.2-inch large backlit screen

The love for Mom and Dad is clearly visible.

## Technical Parameters

Measuring mode:	Step-down measurement (oscillometric method)
Measurement accuracy:	Static pressure $\pm 3$ mmHg (0.4 kPa) Pulse: 5%
Measuring range:	Blood pressure (0-280 mm Hg/0-37.3kpa); Pulse (40-199 beats/min)
Operating environment:	Temperature (5°C to 40°C ); humidity (15% to 85%RH)
Storage environment:	Temperature (-20°C to 55°C ); humidity ( $\leq 93\%$ RH)
Storage capacity:	90 sets of data each for two users
Main unit size:	126.5*126.5*52mm
Single unit weight:	407g (including battery)

### Shenzhen Lepu Intelligent Medical Device Co., Ltd.

Address: Floor 3, BLD 9, BaiWangxin High-Tech Industrial Park,  
Songbai Road, Xili Street, Nanshan District, Shenzhen, China.

E-mail: paul.han@lepu-medical.com

Web: www.lepucare.com

Tel: 400-830-9392



# Fingertip Pulse Oximeter

## INTRODUCTION:

Lepu Pulse Oximeter delivers easy, fast and reliable measurements for SpO<sub>2</sub>, PR and PI with its compact and light body.

It can be widely used in families, hospitals, community healthcare centers and health clubs etc.



Oxygen Saturation

Pulse Rate



PCTG material



# Fingertip Pulse Oximeter

## MAIN FEATURES:

- OLED color display, easy to read
- SpO<sub>2</sub>, PR, PI and bar graph can be displayed on the screen
- Automatically power off within 8 seconds when no finger detected
- Low power indicator
- 4 Kinds of color available (Blue, Green, Cyan and Yellow)



Blue  
LOX100A

Green  
LOX100B

Cyan  
LOX100C

Yellow  
LOX100D

## SPECIFICATIONS:

Oxygen Saturation	Measuring Range	70%~99%
	Accuracy	80%~99%, ±2% 70%~79%, ±3% No requirement for below 70%
	Resolution	1%
Pulse Rate	Measuring Range	30bpm~240bpm
	Accuracy	30bpm~240bpm, ±2bpm or ±2% (select larger)
	Resolution	1bpm
Blood Perfusion Index	Measuring Range	0.3%~20%
Battery Model		2 AAA Batteries
Power Consumption		< 30mA
Battery Life		Continuous use for 25 hours

### Shenzhen Lepu Intelligent Medical Equipment Co.,Ltd.

Floor 3, BLD 9, BaiWangxin High-Tech Industrial Park, SongbaiRoad, Xili Street, Nanshan District, Shenzhen, China



Company name: Lepu Medical (Europe) Cooperatief U.A.

Add: AbeLenstra Boulevard 36, 8448

JB, Heerenveen, The Netherlands

Http: //www.lepuhomecare.com

E-mail: paul.han@lepu-medical.com



**+86-755-86952278**

LEPU 乐普

# LEPU LFR30B

Infrared forehead thermometer



CE 2797 FDA



Non-Contact



One Second Measurement



IP22 Protection



3 Colors Backlights



Indication for fever



99 Sets reading Data



- ① Memory button
- ② Display
- ③ Sensing probe
- ④ Mode switch button
- ⑤ On/Off button
- ⑥ Bottom shell

## Product Information

- Brand: LEPU
- Model: LFR30B
- Product name: LEPU Infrared Forehead Thermometer
- Size: 168.5(L) mm x 36(W)mm x 48(H) mm
- Weight: 95g (including battery)
- Measurement range: 32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
- Laboratory Accuracy : Between 34°C and 43°C, requires  $\pm 0.3^{\circ}\text{C}$  Not within this range, requires  $\pm 0.4^{\circ}\text{C}$
- Operating temperature: 5°C ~ 40 °C

### Shenzhen lepu intelligent medical device co. LTD

📍 Add: Floor 3, BLD 9, BaiWangxin High-Tech Industrial Park,  
Songbai Road, Xili Street, Nanshan District, Shenzhen, China.

✉ E-mail: paul.han@lepu-medical.com

🌐 Web: www.lepucare.com

☎ Tel: 400-830-9392





LEPU 乐普

# LEPU LFR50

Infrared Forehead Thermometer



CE<sub>2797</sub> FDA



One Second Measurement



Non-Contact



IP22 Protection



3 Colors Backlights



Indication for fever



99 Sets reading Data



## Product Information

- Brand: LEPU
- Model: LFR50
- Product name: LEPU Infrared Forehead Thermometer
- Size: 132\*42\*173.5mm
- Weight: 192g (including battery)
- Measurement range: 32.0 °C ~ 43.0 °C
- Laboratory Accuracy : Between 34°C and 43°C requires  $\pm 0.3^{\circ}\text{C}$  Not within this range, requires  $\pm 0.4^{\circ}\text{C}$
- Operating temperature: 5°C ~ 40 °C

### Shenzhen lepu intelligent medical device co. LTD

📍 Address: North side of floor 3, BLD 9 BaiWangxin High-Tech Industrial Park Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

✉ E-mail: paul.han@lepu-medical.com

🌐 Website: www.lepucare.com

☎ Tel: 400-830-9392



LEPU 乐普

# LEPU LFR60

Infrared Forehead Thermometer



CE 2797

FDA



Available in low-temperature environment



One Second Measurement



IP22 Protection




3 Colors Backlights



Indication for fever



99 Sets reading Data



1.Sensor probe



2.Display



3.Mode button



4.Memory button



5.On/Offbutton



6.Bottom shell

## Product Information

- Brand: LEPU
- Model: LFR60
- Product name: LEPU Infrared Forehead Thermometer
- Size: 129\*41.5\*171mm
- Weight: 188g (including battery)
- Measurement range: 32.0 °C ~ 43.0 °C
- Accuracy:  $\pm 0.3$  °C (32 °C ~ 43 °C)  $\pm 0.4$  °C (Not within 32 °C ~ 42 °C)
- Operating temperature: 5 °C ~ 40 °C

### Shenzhen lepu intelligent medical device co. LTD

📍 Address: North side of floor 3, BLD 9 BaiWangxin High-Tech Industrial Park Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

✉ E-mail: paul.han@lepu-medical.com

🌐 Website: www.lepucare.com

☎ Tel: 400-830-9392



# LMT 10 Digital Thermometer

LMT 10 digital thermometer has a pure white light body. The light weight and small size make the product easy to use and carry. The digital thermometer can measure the human body temperature quickly and accurately with the buzzing reminder. In addition to high measurement accuracy, LMT 10 is also a device with IP22 waterproof level.

Lightweight

Reminder

Fast & Accurate



# Performance Parameters

## Standard Type: Medical Standard



Division Value:  
0.1° C / 0.1° F



Measurement Range:  
32.0° C to 42.9° C (89.6° F to 109.2° F)



Product Size:  
130.5(length) mm X 20.5  
(width) mm X 11.4(height) mm



Accuracy:  
(35°C~42.0°C) ±0.1°C;  
(32.0~34.9° C and 42.1~42.9° C) ± 0.2° C



Buzzing Reminder:  
Yes



Product Weight:  
12g (Including Battery)



Can Be Switched  
Between ° C and ° F:  
Yes



Water Proof:  
Yes



Battery Lifetime:  
Approx. 10000 Measurements  
(Using a New Battery)



Battery:  
LR41 (1.5V) / SR41 (1.55V)



Operating Conditions:  
5-40° C



Automatic Shutdown Time:  
8 Minutes 40 Seconds ± 1 Minute



Storage / Transport Conditions:  
-20 ~ +55° C



Normal Body Temperature Range:  
Armpit temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278 Fax: +86 0755-86952278 ✉ E-mail: paul.han@lepu-medical.com

# LMT 11A Digital Thermometer

LMT 11A digital thermometer has the light body with blue and white. Its light weight makes it easy to use and carry. The intelligent reminder function makes it convenient for people to use it, and LMT 11A can be also switched the temperature between Fahrenheit and Celsius.

Lightweight

Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard



Division Value:  
0.1° C / 0.1° F



Measurement Range:  
32.0° C to 42.9° C (89.6° F to 109.2° F)



Product Size:  
131.8(Length) mm X 21.1  
(Width) mm X 11.9(Height) mm



Accuracy:  
(35°C~42.0°C) ±0.1°C;  
(32.0~34.9° C and 42.1~42.9° C) ±0.2° C



Buzzing Reminder:  
Yes



Product Weight:  
12g (Including Battery)



Can Be Switched  
Between ° C and ° F:  
Yes



Water Proof:  
Yes



Battery Lifetime:  
Approx. 10000 Measurements  
(Using a New Battery)



Battery:  
LR41 (1.5V) / SR41 (1.55V)



Operating Conditions:  
5-40° C



Automatic Shutdown Time:  
8 Minutes 40 Seconds ± 1 Minute



Storage / Transport Conditions:  
-20 ~ +55° C



Normal Body Temperature Range:  
Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278 Fax: +86 0755-86952278 ✉ E-mail: paul.han@lepu-medical.com



# LMT 11B Digital Thermometer

LMT 11B digital thermometer has the light body with blue and white. Its temperature measuring probe is soft, so the LMT 11B is also suitable for the elderly and babies. The intelligent reminder function makes it convenient for people to use it, and LMT 11B can be also switched the temperature between Fahrenheit and Celsius.

Lightweight

Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard



Division Value:  
0.1 ° C / 0.1 ° F



Measurement Range:  
32.0 ° C to 42.9 ° C (89.6 ° F to 109.2 ° F)



Product Size:  
131.8 (Length) mm X 21.1  
(Width) mm X 11.9 (Height) mm



Accuracy:  
(35 ° C ~ 42.0 ° C) ± 0.1 ° C;  
(32.0 ~ 34.9 ° C and 42.1 ~ 42.9 ° C) ± 0.2 ° C



Buzzing Reminder:  
Yes



Product Weight:  
12g (Including Battery)



Can Be Switched  
Between ° C and ° F:  
Yes



Water Proof:  
Yes



Battery Lifetime:  
Approx. 10000 Measurements  
(Using a New Battery)



Battery:  
LR41 (1.5V) / SR41 (1.55V)



Operating Conditions:  
5-40 ° C



Automatic Shutdown Time:  
8 Minutes 40 Seconds ± 1 Minute



Storage / Transport Conditions:  
-20 ~ +55 ° C



Normal Body Temperature Range:  
Armpit Temperature: 36.0 ° C - 37.3 ° C; Rectal Temperature: 36.5 ° C - 37.6 ° C;  
Oral Temperature: 35.9 ° C - 36 ° C



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278   Fax: +86 0755-86952278   ✉ E-mail: paul.han@lepu-medical.com

# LMT 11C Digital Thermometer

LMT 11C digital thermometer has a lightweight body with gray and white. Buzzing reminder and water - proof function are its two highlights. Whether at home or in the hospital, people can measure body temperature quickly and accurately.

Lightweight















Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

-  **Division Value:**  
0.1° C / 0.1° F
-  **Measurement Range:**  
32.0° C to 42.9° C (89.6° F to 109.2° F)
-  **Product Size:**  
131.8(length) mm X 21.1  
(width) mm X 11.9(height) mm
-  **Accuracy:**  
(35°C~42.0°C) ±0.1°C;  
(32.0~34.9° C and 42.1~42.9° C) ± 0.2° C
-  **Buzzing Reminder:**  
Yes
-  **Product Weight:**  
12g (Including Battery)
-  **Can Be Switched  
Between ° C and ° F:**  
Yes
-  **Water Proof:**  
Yes
-  **Battery Lifetime:**  
Approx. 10000 Measurements  
(Using a New Battery)
-  **Battery:**  
LR41 (1.5V) / SR41 (1.55V)
-  **Operating Conditions:**  
5-40° C
-  **Automatic Shutdown Time:**  
8 Minutes 40 Seconds ± 1 Minute
-  **Storage / Transport Conditions:**  
-20 ~ +55° C
-  **Normal Body Temperature Range:**  
Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278 Fax: +86 0755-86952278 ✉ E-mail: paul.han@lepu-medical.com

# LMT 12A Digital Thermometer

LMT 12A digital thermometer can measure the human body temperature quickly and accurately in 20 seconds. It has the characteristics of convenient reading, short measurement time, high measurement accuracy, buzzer reminder, and no mercury. It is especially suitable for use in hospitals, homes and other environments.

Lightweight












Reminder


Fast & Accurate



# Performance Parameters

## Standard Type: Medical Standard

-  **Division Value:**  
0.1° C / 0.1° F
-  **Product Size:**  
133.5(length) mm X 32.0  
(width) mm X 15.0(height) mm
-  **Buzzing Reminder:**  
Yes
-  **Battery Lifetime:**  
Approx. 10000 Measurements  
(Using a New Battery)
-  **Can Be Switched  
Between ° C and ° F:**  
Yes
-  **Operating Conditions:**  
5-40° C
-  **Storage / Transport Conditions:**  
-20 ~ +55° C
-  **Normal Body Temperature Range:**  
Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C
-  **Measurement Range:**  
32.0° C to 42.9° C (89.6° F to 109.2° F)
-  **Accuracy:**  
(35°C~42.0°C) ±0.1°C;  
(32.0~34.9° C and 42.1~42.9° C) ±0.2° C
-  **Approx. Measuring Time:**  
20 Seconds
-  **Product Weight:**  
22g (Including Battery)
-  **Water Proof:**  
Yes
-  **Battery:**  
CR2032 (3.0V)
-  **Automatic Shutdown:**  
Yes



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278   Fax: +86 0755-86952278   ✉ E-mail: paul.han@lepu-medical.com

# LMT 12B Digital Thermometer

LMT 12B digital thermometer can measure the human body temperature quickly and accurately in 20 seconds. It has the soft temperature probe, the large screen, the buzzing reminder function, high measurement accuracy, and no mercury. It's especially suitable for use in hospitals, homes and other environments.

Lightweight

Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard



Division Value:  
0.1° C / 0.1° F



Measurement Range:  
32.0° C to 42.9° C (89.6° F to 109.2° F)



Product Size:  
133.5(length) mm X 32.0  
(width) mm X 15.0(height) mm



Accuracy:  
(35°C~42.0°C) ±0.1°C;  
(32.0~34.9° C and 42.1~42.9° C) ±0.2° C



Buzzing Reminder:  
Yes



Approx. Measuring Time:  
20 Seconds



Battery Lifetime:  
Approx. 10000 Measurements  
(Using a New Battery)



Product Weight:  
22g (Including Battery)



Can Be Switched  
Between ° C and ° F:  
Yes



Water Proof:  
Yes



Operating Conditions:  
5-40° C



Battery:  
CR2032 (3.0V)



Storage / Transport Conditions:  
-20 ~ +55° C



Automatic Shutdown:  
Yes



Normal Body Temperature Range:  
Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C



Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278 Fax: +86 0755-86952278 ✉ E-mail: paul.han@lepu-medical.com