

PC-66B

Handheld Pulse Oximeter



Features

- 2.2 inch TFT display for SpO2, PI, Plethysmogram, PR
- MoveOxy™ SpO2 technology with artifact removal and anti-motion performance
- Various choice of SpO2 probe for Adult, Pediatric and Neonate particularly
- Inbuilt lithium battery or AA battery for continuous monitoring
- Up to 384 hours SpO2 data storage for review
- Adjustable alarm limit, pitch tone level and key tone



Technical Specifications

SpO2

Technique: Dual-wavelength LED

Measuring range: 35%~100%

Measuring accuracy: $\pm 2\%$ (range from 70% to 100%)

Pulse Rate

Measuring range: 30bpm~250bpm

Measuring accuracy: ± 2 bpm or $\pm 2\%$

Perfusion Index

PI Range: 0.2%~20%

Data Storage

Up to 384 hours

Display

2.2 TFT display

Power Supply

AA disposable battery or Lithium battery

Charging

USB cable with adapter

Dimension

145 mm x 74 mm x 29 mm

Net Weight

210g

Data Upload

USB cable

Operating Environment

Operating Temperature: 5°C ~40°C

Operating Humidity: 30%~80%



Adult Clip Sensor



Adult Soft Sensor



Pediatric Clip Sensor



Pediatric Soft Sensor



Neonate/Infant wrap sensor



Shenzhen Creative Industry Co., Ltd.

Floor 5, BLD 9, Baiwangxin High-Tech Industrial Park, Songbai Road, Xili Street

Nanshan District, 518110 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 755 2643 3514 Fax: +86 755 2643 0390

www.creative-sz.com info@creative-sz.com